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Ⓔ Infrared health lamp

Instruction manual

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
INTRODUCTION

Infrared light stimulates blood circulation and your well-being. Warmth deeply penetrates tissues and muscles. **Using red light does not substitute medical treatment. If you have serious muscle injuries, contusions, persistent cold or other disorders, you should consult your doctor first.**

Red light should not be used for medical treatment, unless stated otherwise by your doctor.

SAFETY ADVICES

1. Read this instruction manual carefully before using the product.
2. If you pass on the device to a third party, please include this instruction manual as well.
3. This device complies with technical norms and safety requirements for electric devices.
4. The manufacturer accepts no responsibility for damages resulting from improper use or non-compliance with this instruction manual.
5. Before using the device check if power supply voltage corresponds with the power supply details on the rating plate.
6. This product is intended only for private, not commercial use.
7. Should damage occur, do not open the device yourself. If the device or its power cord is damaged, consult a qualified specialist.
8. Make sure that the power cord is not rubbed against sharp edges or is not folded. The device must not be used if the power cord is visibly damaged.
9. In order to avoid any hazards, a damaged outer flexible cord of this lamp can be replaced only by the manufacturer, its Service Centre or similarly qualified specialist.
10. Make sure that nobody trips over the power cord.
11. **Do not immerse the device, power cord or power plug in water. Do not operate the device with wet hands. Danger of electric shock!**
12. The device must not be used in the proximity of a bathtub or shower, or a washbasin filled with water.
13. **Should the device fall in water, have it checked by a qualified specialist. Do not touch it with water! Always disconnect the power plug first!**
14. Do not use the device outdoors.
15. This device is not intended for use by persons (including children) with reduced physical, sensory, or mental capabilities, or by those who lack experience and knowledge, unless they are supervised or have been given instructions concerning the use of the device by a person responsible for their safety.
16. Children must be supervised at all times to ensure they do not play with the device.
17. Keep electric devices out of the reach of children.
18. The lamp cannot be used if it has been visibly damaged or should have been dropped.
19. Do not disconnect the power plug by pulling the power cord.
20. Use the device only on a level, non-slip surface. Remove flammable and metallic objects from the area of radiation.
21. Should you suffer from skin swellings or rash, burns, blood circulation disorders, diabetes, tuberculosis, tumours, contusions, neurotic pains or aches of unknown origin, or if you have a pace maker, it is strongly advised not to use the device.
22. Eyes must be protected from infrared radiation.

23. **The lamp must not be covered and has to be maintained dust-free at all times. Danger of fire!**
24. **Do not operate the lamp, or leave it to cool down, behind curtains, in the proximity thereof or other flammable materials. Danger of fire!**
25. It is strongly recommended to keep the minimum distance of 50 cm to  the irradiated surface (also skin).
26. You can use only illuminates of the type PHILIPS infrared light bulb R95E (230V~, 100W). Light bulbs of other types may damage the device.

This device is suitable only for private use and must not be used for commercial purposes.

OPERATION

1. Mount the device following the *Safety Advice*s carefully.
2. Connect the power plug with a corresponding socket.
3. The red light starts glowing. Infrared light heats up very quickly, pre-heating is not necessary.
4. Adjust the angle of the lamp keeping the distance of 50 cm to the irradiated surface.
5. **Do not use the device for more than 15 minutes at a time!**
6. **You can use the device repeatedly, however, you should wait a couple of hours before using the lamp again.**
7. After each use the lamp must cool down completely.
8. Surface of the lamp gets very hot during use. In order to adjust the lamp, touch it only on the rear end of the light head. **Burn hazard!**

CLEANING AND MAINTENANCE

1. Before cleaning, disconnect the power plug from the socket and let the device cool down completely.
2. Clean the lamp only with a dry or slightly dampened cloth.

REPLACING THE LIGHT BULB

1. Should the bulb do not work, it needs to be replaced.
2. Before replacing, disconnect the lamp from the mains supply and let the device cool down completely.
3. Turn the bulb counter clockwise, take it out and dispose of it according to binding regulations.
4. Insert a new bulb of the same type [PHILIPS infrared light bulb R95E (230V~, 100W)].
5. The lamp can be used again.

TECHNICAL DATA

Red light: Max. 100W
Rated voltage: 230V~ 50Hz
Light bulb: PHILIPS infrared light bulb R95E (230V~, 100W)

Environmental protection

Electric and electronic devices should not be disposed of with household garbage. If you want to dispose of this device, contact local authorities about the correct method of waste disposal.

When the operational life of the device has come to an end, disable the device by disconnecting the power plug from the socket and cut the power cord in two.

This device is suitable only for private use and must not be used for commercial purposes.

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Warranty

(Infrared health lamp / EAN 430561503062 3)

We guarantee the perfect quality and functionality of our technical appliances pursuant to our warranty conditions. The warranty period is 24 months and begins on the date of purchase. Please keep your receipt as proof of purchase. If you should have reason to complain, despite our extensive quality control procedures, please call the free service number stated below as the first step. The line is staffed by qualified employees who will make a note of your problem, give you a transaction number and arrange for replacement goods if necessary. We shall send you a label for the return consignment by e-mail, fax or post, which you can use to send the product to our service partner free of charge. We regret that any appliances which are sent in without such a transaction number cannot be accepted. The warranty does not apply to damage caused by improper use of the appliance, wearing parts or consumables. If possible, please use the original packaging to return the goods to avoid damage in transit.

This appliance has undergone strict quality control to ensure the perfect quality of the materials, workmanship and function. If, however, a fault should occur which is not attributable to improper use of the product (e.g. caused by water or dropping the appliance, incorrect operation etc.), we undertake to remedy the damage under the warranty (by replacing or repairing faulty parts). Within the scope of its warranty obligations, Rossmann shall remedy only damage to the product directly covered by the warranty, in the form of replacement delivery or repair.

The above warranty does not cover further claims for compensation of any kind, in particular claims owing to breach of obligations or tort, claims for reimbursement of outlay incurred to no avail, whether these claims are filed against Rossmann itself or against any agents employed in performance of Rossmann's obligations; the same applies to compensation for consequential damage. No further claims whatsoever shall be created within the course of the warranty. The above restrictions do not apply to the statutory rights to which the purchaser (consumer) is entitled owing to the purchase of consumer goods, to damage arising from injury to life, limb or health, to damage which is attributable to intent or gross negligence or to claims pursuant to the Product Liability Act. The warranty ceases to apply in case of improper use, unauthorised intervention or any technical modification of the appliance by the consumer. Any performance under the warranty does not constitute extension or renewal of the warranty period of 24 months. Replacement batteries are explicitly exempt from the warranty.

An exact description of your complaint will enable us to repair it more quickly. Even if the warranty has expired, you still have the option of calling the Ideenwelt service partner hotline. However, you will be charged for any repairs required after expiry of the warranty period, as these are no longer covered by the warranty.

Jet GmbH Hotline: 01805 - 2288904

Please complete in block capitals and return with the product:

Name _____

Address _____

Postcode, city _____

Telephone (daytime) _____

Date, signature _____

Description of the article _____

Description of the damage _____

If the warranty does not apply
(please tick as appropriate)

- ☐ Please return the product unrepaid. I agree to assume the cost of transport.
- ☐ Please repair the product and send me the invoice.

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